Product datasheet

MON9826



Mouse anti-Haemoglobin, clone EBS-O-048 (Monoclonal) Clone no. EBS-O-048

Product name	Mouse anti-Haemoglobin, clone EBS-O-048 (Monoclonal)
Host	Mouse
Applications	ELISA, FC, IHC-fr, IF
Species reactivity	human
Conjugate	-
Immunogen	purified human haemoglobin
lastime	
lsotype	lgG3
Clonality	Monoclonal
Clone number	EBS-O-048
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON9826



Mouse anti-Haemoglobin, clone EBS-O-048 (Monoclonal) Clone no. EBS-O-048

MONOSAN

Additional info

EBS-O-148 can be used in ELISA, frozen sections, FACS and Western blot, however no method is known to date to retrieve the antigen in paraffin sections. Anti-hemoglobin antibodies can be used in diagnosis of anemias and as tumor marker in stool.

1.	-
2	-
3.	-
4.	-
5.	-
	4.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com