Product datasheet

MON9831



Rat anti-Laminin y, clone A5 (Monoclonal) Clone no. A5

Product name	Rat anti-Laminin y, clone A5 (Monoclonal)
Host	Rat
Applications	FC, IHC-fr, IHC-P, IF
Species reactivity	human, monkey, mouse, pig, rat, rhesus
Conjugate	-
Immunogen	laminin preparation from murine EHS cells
lsotype	lgG2a-K
Clonality	Monoclonal
Clone number	A5
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON9831



Rat anti-Laminin y, clone A5 (Monoclonal) Clone no. A5

MONOSAN

Additional info

Laminins are large hetero-trimeric, non-collagenous glycoproteins found in basement membranes and composed of α , β , and γ chains. A5 reacts specifically with γ chain 1. Alterations of basement membrane integrity, from local discontinuities up to complete loss, are described in many types of human and animal epithelial neoplasms.

References	1. 2	Ljubimov JY et. al. Cancer Res. 61(14): 5601-5610 (2001) Ljubimov AV et. al. Int J Cancer 50: 562-566 (1992)
	3.	Ljubimov AV et. al. Lab Invest 72: 461-473 (1995)
	4.	Folkvord et. al. J Histochem Cytochem 37: 105-113 (1989)
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021