Product datasheet MON9837



## Mouse anti-Nucleolar Antigen, clone AE-3 (Monoclonal)

Clone no. AE-3 MONOSAN

Product name Mouse anti-Nucleolar Antigen, clone AE-3 (Monoclonal)

**Host** Mouse

**Applications** IHC-fr, IHC-P, IF

Species reactivity human

Conjugate -

Immunogen nuclei of myeloid leukemia biopsy cells

Isotype IgG1-K

**Clonality** Monoclonal

Clone number AE-3

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9837



Mouse anti-Nucleolar Antigen, clone AE-3 (Monoclonal)

Clone no. AE-3 MONOSAN

Additional info

AE-3 recognizes a nucleolar antigen expressed in nucleoli of human cells. AE-3 can also be used as a marker of the nucleolus in subcellular fractions.

References 1. Epstein, AL et al. Progress on non-histone protein res. 1:117-137 (1985)

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES