Product datasheet MON9841



## Mouse anti-TGF-alpha, clone 1E8-G6 (Monoclonal)

Clone no. 1E8-G6 MONOSAN

Product name Mouse anti-TGF-alpha, clone 1E8-G6 (Monoclonal)

**Host** Mouse

**Applications** ELISA, FC, IHC-fr, IHC-P, WB

Species reactivity human, rabbit, zebrafish

Conjugate -

Immunogen synthetic TGF-alpha 34-43 decapeptide linked carrier

Isotype IgM-K

**Clonality** Monoclonal

Clone number 1E8-G6

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9841



Mouse anti-TGF-alpha, clone 1E8-G6 (Monoclonal)

Clone no. 1E8-G6 MONOSAN

## Additional info

1E8-G6 reacts with TGF-alpha and shows no cross-reaction with EGF and the neuropeptide synenkephalin. 1E8-G6 is completely blocked by the peptide used for immunization. TGFalpha (aa50) is a growth factor with 33% homology to EGF, binds to EGFR, activates tyrosine phosphorylation of the receptor, and stimulates cell proliferation. It plays a role in tumor initiation by inducing the reversible transformed phenotype. In sections both cytoplasma and cell membranes are stained

References 1. Bebók, Zs, et al. Breast Cancer Res.Treatm. 29: 229-235 (1994)

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES