## Product datasheet

MON9847



Mouse anti-IgG, clone EBS-huG (Monoclonal)Clone no.EBS-huG

Product name	Mouse anti-IgG, clone EBS-huG (Monoclonal)	
Host	Mouse	
Applications	ELISA, FC, IHC-fr	
Species reactivity	human	
Conjugate	-	
Immunogen	purified human immunoglobulins	
lsotype	lgG2a-K	
Clonality	Monoclonal	
Clone number	EBS-huG	
Size	100 ug	
Concentration	100 ug/ml	
Format	-	
Storage buffer	PBS with 0.02% sodium azide	
Storage until expiry date	2-8°C	

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Product datasheet

MON9847



Mouse anti-IgG, clone EBS-huG (Monoclonal) Clone no. EBS-huG

Additional info

EBS-huG reacts specifically with human IgG and not with other types of immunoglobulins.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com