

Mouse anti-Keratin 18, clone C-04 (Monoclonal)

Clone no. C-04

MONOSAN

Product name	Mouse anti-Keratin 18, clone C-04 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P, IF, WB
Species reactivity	cat, cow, dog, ferret, goat, hamster, horse, human, marmoset, mouse, pig, rat, :
Conjugate	-
Immunogen	cytoskeleton preparation of human A431 carcinoma cell line
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	C-04
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Keratin 18, clone C-04 (Monoclonal)

Clone no. C-04

MONOSAN

Additional info

C-04 Reacts with keratin 18 (45kDa) present in a wide variety of simple epithelia. It does not react with stratified squamous epithelia. It reacts with epithelial tumors of the gastrointestinal tract, lung, breast, pancreas, ovary, and thyroid.

References

1. Bártek et al. J. of Pathol. 164: 215-24 (1991)
2. Bártková et al. Neoplasma 38(4): 439-46 (1991)
3. Lauerová et al. Hybridoma 7: 495-504 (1988)
4. Taylor –Papadimitriou et.al. J. Cell Sci. 94: 403-13 (1989)
5. Kovařík et al. J. Tumor Marker Oncol. 5: 219 (1990)

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES