Product datasheet MON9861



Mouse anti-Keratin 4,5,6,8,10,13,18, clone C-11 (Monoclonal)

Clone no. C-11 MONOSAN

Product name Mouse anti-Keratin 4,5,6,8,10,13,18, clone C-11 (Monoclonal)

Host Mouse

Applications FC, IHC-fr, IHC-P, IF, WB

Species reactivity cat, chicken, cow, dog, ferret, goat, guinea pig, human, kangaroo, marmoset, m

Conjugate -

Immunogen cytoskeleton preparation of human A431 carcinoma cell line

Isotype IgG1-K

Clonality Monoclonal

Clone number C-11

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9861



Mouse anti-Keratin 4,5,6,8,10,13,18, clone C-11 (Monoclonal)

Clone no. C-11 MONOSAN

Additional info

C11 reacts with keratins: 4,5,6,8,10,13 and 18. This is a broad-spectrum antibody, which has been reported to differentiate epithelial tumors from non-epithelial tumors. Many studies have shown the usefulness of keratins as markers in cancer research and tumor diagnosis.

References 1. Vojtěšek, B. et al. Folia Biol., 35(6): 373-382, (1989)

2 Vojtěšek, B. et al. Neoplasma 37(3): 333-342 (1990)

3. Kovařík J. et al. Int. J. Cancer, Suppl. 3: 50-55, (1988)

4. Kovařík J. et al. J. Tumor Marker Oncol. 5: 219 (1990)

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES