

Mouse anti-MHC I Aw32+A25, clone CATA-1 (Monoclonal)

Clone no. CATA-1

MONOSAN

---

Product name	Mouse anti-MHC I Aw32+A25, clone CATA-1 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P
Species reactivity	human
Conjugate	-
Immunogen	human PBL of phenotype A1, Aw32, B7, B37, CW-, DR2, DRw10
Isotype	IgG2a-K
Clonality	Monoclonal
Clone number	CATA-1
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MHC I Aw32+A25, clone CATA-1 (Monoclonal)

Clone no. CATA-1

MONOSAN

**Additional info**

CATA-1 reacts with cells bearing HLA-A25 or HLA-Aw32 antigens. In addition, a reaction was observed with a cell of phenotype A2, Aw31; B17, Bw49.

**References**

1. Vilella, R. et al. Human Immunol. 6: 53-62 (1983)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES