

Mouse anti-MHC II DRA, clone EBS-O-111 (Monoclonal)

Clone no. EBS-O-111

MONOSAN

Product name	Mouse anti-MHC II DRA, clone EBS-O-111 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P
Species reactivity	human
Conjugate	-
Immunogen	RAJI cells
Isotype	IgG2a-K
Clonality	Monoclonal
Clone number	EBS-O-111
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MHC II DRA, clone EBS-O-111 (Monoclonal)

Clone no. EBS-O-111

MONOSAN

Additional info

MHC class II molecules are encoded by polymorphic MHC genes and consist of a non-covalent complex of an α and β chain. Helper T lymphocytes bind antigenic peptides presented by MHC class II molecules. MHC class II molecules bind 13-18 amino acid antigenic peptides. Accumulating in endosomal/lysosomal compartments and on the surface of B-cells, HLA-DM and -DO molecules regulate binding of exogenous peptides to class II molecules (HLA-DR) by sustaining a conformation that favors peptide exchange. The differential structural properties of MHC class I and class II molecules account for their respective roles in activating different populations of T lymphocytes

References

1. Thompson CJ et al., Human Immunol 6: 133-150 (1983)
2. Rask L, et al., Autoimmunity 8(3):237-244 (1991)
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES