

Mouse anti-Cytomegalovirus p65, clone CMV-227 (Monoclonal)

Clone no. CMV-227

MONOSAN

Product name	Mouse anti-Cytomegalovirus p65, clone CMV-227 (Monoclonal)
Host	Mouse
Applications	ELISA, FC, IHC-fr, IF
Species reactivity	CMV
Conjugate	-
Immunogen	Recombinant CMV-p65 protein
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	CMV-227
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cytomegalovirus p65, clone CMV-227 (Monoclonal)

Clone no. CMV-227

MONOSAN

Additional info

Ten to forty percent of the patients with acquired immunodeficiency syndrome (AIDS) develop cytomegalovirus (CMV) infections. In some patients with AIDS, CMV is detected in the bronchoalveolar lavage fluid (BALF), urine, and other specimens, even when there are no symptoms of CMV disease. An indicator of active CMV infection is needed to facilitate the diagnosis of CMV disease in patients with AIDS or HIV infection. CMV p65 antigen was detected in the leukocytes of both the peripheral blood and BALF during the early phase of CMV disease.

References

1. Hanabusa H, et al, Kansenshogaku Zasshi. 68(9): 1105-12 (1994)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES