

Mouse anti-EBV-LMP1, clone EBS-I-025 (Monoclonal)

Clone no. EBS-I-025

MONOSAN

Product name	Mouse anti-EBV-LMP1, clone EBS-I-025 (Monoclonal)
Host	Mouse
Applications	ELISA, IF, IP, IHC-P, WB
Species reactivity	EBV
Conjugate	-
Immunogen	LMP1 b-galactosidase fusion protein (AA 198-486)
Isotype	IgG1
Clonality	Monoclonal
Clone number	EBS-I-025
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-EBV-LMP1, clone EBS-I-025 (Monoclonal)

Clone no. EBS-I-025

MONOSAN

Additional info

EBS-I-025 reacts with a conserved repeat domain on LMP1 (AA 268-286), present in all EBV isolates. LMP1 is expressed in most viral latency stages in vitro and in vivo

References

1. Meij, P. et al. J Infect Dis. 179: 1108-15 (1999)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES