Product datasheet MON9891



Mouse anti-HPV-16 E6, clone HPV-13E2 (Monoclonal)

Clone no. HPV-13E2 MONOSAN

Product name Mouse anti-HPV-16 E6, clone HPV-13E2 (Monoclonal)

Host Mouse

Applications ELISA, IHC-fr

Species reactivity human

Conjugate -

Immunogen HPV-16 E6 synthetic peptide

Isotype IgG1-K

Clonality Monoclonal

Clone number HPV-13E2

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9891



Mouse anti-HPV-16 E6, clone HPV-13E2 (Monoclonal)

Clone no. HPV-13E2 MONOSAN

Additional info

Monoclonal antibody HPV-13E2 reacts with HPV-16 E6 peptide (MHQKRTAMFQDPPQERPRKLPQLC). Infection with specific types of HPV has been associated with an increased risk of developing cervical neoplasia. HPV types 6 and 11 have been associated with relatively benign diseases such as genital warts but types 16 and 18 are strongly associated with cervical, vaginal, and vulvar malignancies

References 1. Beldermann, F. Virusprotein E6 des humanen Papillomavirus HPV-16, Thesis, L

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES