

Mouse anti-BCL 10, clone EBS-T-001 (Monoclonal)

Clone no. EBS-T-001

MONOSAN

Product name	Mouse anti-BCL 10, clone EBS-T-001 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P, WB
Species reactivity	human
Conjugate	-
Immunogen	human BCL10 recombinant protein (aa122-168)
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	EBS-T-001
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-BCL 10, clone EBS-T-001 (Monoclonal)

Clone no. EBS-T-001

MONOSAN

Additional info

EBS-T-001 reacts with BCL10. Having an N-terminal caspase recruitment domain (CARD), BCL10 can induce apoptosis and activate NF- κ B. It is found on subpopulations of normal B and T cells, and is associated with MALT1, a paracaspase that, like BCL10, can be found translocated in MALT lymphoma. In such cases either BCL10 or MALT1 or both are highly expressed, depending on the site of translocation. MALT lymphomas lacking this translocation exhibit much lower levels of expression. BCL10 has been shown to be functionally conserved all the way back to zebrafish

References

1. Ye H, et al, Am J Pathol 157: 1147-54 (2000)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES