Product datasheet MON9898



## Mouse anti-BLA36, clone EBS-T-002 (Monoclonal)

Clone no. EBS-T-002 MONOSAN

Product name Mouse anti-BLA36, clone EBS-T-002 (Monoclonal)

**Host** Mouse

**Applications** IHC-fr

Species reactivity human

Conjugate -

Immunogen human BLA36 recombinant protein

Isotype IgG1-K

**Clonality** Monoclonal

Clone number EBS-T-002

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9898



Mouse anti-BLA36, clone EBS-T-002 (Monoclonal)

Clone no. EBS-T-002 MONOSAN

## Additional info

BLA36 is a developmentally regulated 36 kDa antigen expressed on the plasma membrane of B lymphocytes, Reed-Sternberg, and mononuclear Hodgkin's cells. It also gives strong staining of B cell lymphomas including follicular center cell lymphomas (large and small cell types), mantle zone lymphomas, and immunoblastic lymphomas.

References 1. -

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES