Product datasheet MON9931



Mouse anti-CD20, clone 93-1B3 (Monoclonal)

Clone no. 93-1B3 MONOSAN

Product name Mouse anti-CD20, clone 93-1B3 (Monoclonal)

Host Mouse

Applications FC, IHC-fr, FUNC

Species reactivity human

Conjugate -

Immunogen stimulated human leucocytes

Isotype IgG1-K

Clonality Monoclonal

Clone number 93-1B3

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9931



Mouse anti-CD20, clone 93-1B3 (Monoclonal)

Clone no. 93-1B3 MONOSAN

Additional info

93-1B3 binds with CD20 which is a 30/33 kDa non-glycosylated transmembrane phosphoprotein with three extensive hydrophobic regions. CD20 is involved in regulation of Bcell activation. It is expressed on the surface of all B-cells beginning at the pro-B phase (CD45R+, CD117+) and progressively increasing in concentration until maturity. Plasma cells are negative. CD20 is retained on many B-cell malignancies. CD20 positive cells are also sometimes found in cases of Hodgkin's disease, myeloma, and thymoma. 93-1B3 has been clustered at the IIIrd and Vth HLDA Workshops.

References 1. Cobbold, S. et al., in leucocyte typing III, Oxford University Press (1987)

2 Schlossman SF et al. Leukocyte Typing V Oxford University Press: 342-352 (199

3. -

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES