

Mouse anti-CD23, clone EBS-CD-016 (Monoclonal)

Clone no. EBS- CD-016

MONOSAN

Product name	Mouse anti-CD23, clone EBS-CD-016 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IF
Species reactivity	human
Conjugate	-
Immunogen	human PBL
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	EBS- CD-016
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD23, clone EBS-CD-016 (Monoclonal)

Clone no. EBS- CD-016

MONOSAN

Additional info

EBS-CD-016 reacts with CD23, which is Fc epsilon Receptor II, which plays a role in the differentiation of B-cells and is B-cell specific.

References

1. Johnson GD, et al, Clin Exp Immuno. 64(1): 205-13 (1986)
2. Bonnefoy JY, et al, J Immunol. 138: 2970 (1987)
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES