Product datasheet MON9974



Mouse anti-CD27, clone 203-6 (Monoclonal)

Clone no. 203-6 MONOSAN

Product name Mouse anti-CD27, clone 203-6 (Monoclonal)

Host Mouse

Applications FC, IHC-fr, IF, WB

Species reactivity human

Conjugate -

Immunogen Sézary cells

Isotype IgG3-K

Clonality Monoclonal

Clone number 203-6

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9974



Mouse anti-CD27, clone 203-6 (Monoclonal)

Clone no. 203-6 MONOSAN

Additional info

203-6 Reacts with human CD27, a disulphide-linked 120 kDa dimer. CD27 is a lymphocyte-specific member of the TNF-R/NGF-R superfamily, and is expressed on a subset of human thymocytes and on the majority of mature T-lymphocytes. CD27 is highly expressed on activated T and B-cells. 203- 6 was clustered at the VIth WLDA.

References 1. Kishimoto T. et al., eds. Leukocyte Typing VI, p509-514, Garland Publishing, Inc,

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES