Product datasheet MON9980



TCC, Human, mAb, clone aE11 (Monoclonal)

Clone no. aE11 MONOSAN

Product name TCC, Human, mAb, clone aE11 (Monoclonal)

Host Mouse

Applications IHC-fr,FC,FUNC,ELISA,IF,IHC-P

Species reactivity human, horse, pig

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype IgG2a

Clonality Monoclonal

Clone number aE11

Size 1 ml

Concentration 100 ug/ ml

Format -

Storage buffer PBS with 0.1% BSA and 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES



TCC, Human, mAb, clone aE11 (Monoclonal)

Clone no. aE11 MONOSAN

Additional info

Monoclonal antibody aE11 reacts with a C9 neoantigen of the terminal complement complex (TCC). The three distinct activation pathways of complement converge with the formation of a C5 convertase. The cleavage of C5 by this convertase initiates the lytic or terminal pathway. In contrast to the activation pathways, which require enzymatic cleavage for activation, the terminal pathway relies on conformational changes induced by binding. Binding of C6 facilitates binding of C7 which alters the conformation of the complex. After binding of C8, a variable number of C9 molecules associate with the C5b678 complex, which is also termed- the terminal complement complex (TCC). The formation of TCC causes lysis of cells or can trigger a variety of cellular metabolic pathways resulting in the synthesis and release of inflammatory mediators. The TCC contains neoantigens that are absent from the individual native components.- C9 neoantigens are present both in the membrane-bound (MAC) and the fluid-phase (SC5b-9) complex. TCC is present in normal human plasma and increased in patients with complement activation.

References

- 1. Mollnes; T et al. Scand | Immunol 1985; 22: 197
- Mollnes, T et al Scand J Immunol 1985, 22: 183
- 3. Pettersen; H et al. Scand | Immunol 1987; 25: 567
- 4. Berstad A et al. Gut 1997; 40:196
- 5. Stewart M et al. BJH 1997; 96:451

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES