Product datasheet MON9990



## Rabbit anti-Uroplakin II (Polyclonal)

Clone no. - MONOSAN

Product name Rabbit anti-Uroplakin II (Polyclonal)

**Host** Rabbit

Applications IHC-P (1:50-1:100)

Species reactivity human

Conjugate -

**Immunogen** Peptide-KLH conjugate

Isotype -

**Clonality** Polyclonal

Clone number -

Size 1 ml

**Concentration** lot specific

Format -

Storage buffer Tissue culture supernatant with sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9990



Rabbit anti-Uroplakin II (Polyclonal)

Clone no. - MONOSAN

## Additional info

Uroplakins (UPs) are a family of transmembrane proteins (UPs Ia, Ib, II and III) that are specific differentiation products of urothelial cells. Uroplakins are markers of terminally differentiated urothelium. Uroplakin II (UPII) is a newly described sensitive marker for urothelial carcinoma (UC). The expression profile of UPII in different types of UC and its utility in the diagnostic setting are needed.

**References** 1. Moll et al. Am J Pathol 1995;147:1383-1397

2 Olsburgh et al. J Pathol 2003;199:41-49

3. Parket et al. 2003;27:1-10

4. Ohtsuka et al. 2006;97:1322-1326

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES