

Mouse anti-Villin, clone CWWB1 (monoclonal)

Clone no. CWWB1

MONXtra

Product name	Mouse anti-Villin, clone CWWB1 (monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr, WB
Species reactivity	human
Conjugate	-
Immunogen	Prokaryotic recombinant protein corresponding to the C-terminal "headpiece" region of the human villin molecule.
Isotype	IgG1, kappa
Clonality	Monoclonal
Clone number	CWWB1
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Villin, clone CWWB1 (monoclonal)

Clone no. CWWB1

MONXtra

Additional info

Analyte Specific Reagent. Analytical and performance characteristics are not established.

References

1. Toyoshima K et al. The Journal of Histochemistry and Cytochemistry. 46 (11): 1.
2. Friederich E et al. BioEssays. 12 (9): 403–408 (1990)
3. Janmey PA and Matsudaira P T. The Journal of Biological Chemistry. 263 (32): 1
4. West AB et al. Gastroenterology. 94: 343–352 (1988)
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES