

Mouse anti-Alpha-Sarcoglycan (Adhalin), clone Ad1/20A6 (monoclonal)

Clone no. Ad1/20A6

MONXtra

Product name	Mouse anti-Alpha-Sarcoglycan (Adhalin), clone Ad1/20A6 (monoclonal)
Host	Mouse
Applications	IHC-fr
Species reactivity	human, mouse, rat, rabbit, hamster, pig
Conjugate	-
Immunogen	Fusion protein containing amino acids 217 to 289 of the rabbit adhalin sequence.
Isotype	IgG1
Clonality	Monoclonal
Clone number	Ad1/20A6
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Alpha-Sarcoglycan (Adhalin), clone Ad1/20A6 (monoclonal)

Clone no. Ad1/20A6

MONXtra

Additional info

Human alpha-sarcoglycan, also known as adhalin. Also crossreacts strongly with alpha-sarcoglycan in sections of muscle from mouse, rat, rabbit, hamster and pig. Does not react with chicken muscle.

References

1. Marafioti T et al. American Journal of Pathology. 162 (3): 861–871 (2003)
2. Hess J et al. Molecular and Cellular Biology. 21 (5): 1531–1539 (2001)
3. Re D et al. Cancer Research. 61 (5): 2080–2084 (2001)
4. Luo Y and Roeder R G. Molecular and Cellular Biology. 15 (8): 4115–4124 (1995)
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES