

Mouse anti-Beta-Sarcoglycan, clone betaSarc1/5B1 (monoclonal)

Clone no. bSarc1/5B1

MONXtra

Product name	Mouse anti-Beta-Sarcoglycan, clone betaSarc1/5B1 (monoclonal)
Host	Mouse
Applications	IHC-fr
Species reactivity	human
Conjugate	-
Immunogen	Fusion protein RBSG-NT of the human beta-sarcoglycan sequence.
Isotype	IgG1
Clonality	Monoclonal
Clone number	bSarc1/5B1
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MONX10814

MONOSAN[®]

Mouse anti-Beta-Sarcoglycan, clone betaSarc1/5B1 (monoclonal)

Clone no. bSarc1/5B1

MONXtra

Additional info

Analyte Specific Reagent. Analytical and performance characteristics are not established

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES