

Mouse anti-Helicobacter pylori, Clone ULC3R (monoclonal)

Clone no. ULC3R

MONXtra

Product name	Mouse anti-Helicobacter pylori, Clone ULC3R (monoclonal)
Host	Mouse
Applications	IHC-P
Species reactivity	human, bacteria
Conjugate	-
Immunogen	Prokaryotic recombinant protein corresponding to 151 amino acids of the HopE molecule.
Isotype	IgG1
Clonality	Monoclonal
Clone number	ULC3R
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with 15mM Sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Helicobacter pylori, Clone ULC3R (monoclonal)

Clone no. ULC3R

MONXtra

Additional info

Tumor cells of epithelial, lymphoid, glial and mesenchymal origin are reported to be negative. This clone is well described in the literature. It is indicated to label an intracytoplasmic antigen in the majority of melanomas and other tumors demonstrating melanoma/melanocytic differentiation.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES