

Mouse anti-B Cell Specific Octamer Binding Protein-1 (BOB-1), clone TG14 (monoclonal)

Clone no. TG14

MONXtra

Product name	Mouse anti-B Cell Specific Octamer Binding Protein-1 (BOB-1), clone TG14 (monoclonal)
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	A prokaryotic recombinant protein corresponding to 116 amino acids of the C-terminal region of the B cell specific octamer binding protein-1 molecule.
Isotype	IgG2b
Clonality	Monoclonal
Clone number	TG14
Size	1ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-B Cell Specific Octamer Binding Protein-1 (BOB-1), clone TG14 (monoclonal)

Clone no. TG14

MONXtra

Additional info

Analyte Specific Reagent. Analytical and performance characteristics are not established.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES