## Product datasheet

MONX11151



MONXtra

Mouse anti-Cytokeratin 8, clone TS1 (monoclonal) Clone no. TS1

Product name	Mouse anti-Cytokeratin 8, clone TS1 (monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr
Species reactivity	human
Conjugate	-
Immunogen	Purified tumor-derived cytokeratins.
lsotype	lgG1
Clonality	Monoclonal
Clone number	TS1
Size	1ml
Concentration	n/a
Format	-
Storage buffer	Tissue culture supernatant with 15mM sodium azide
Storage until expiry date	2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Product datasheet

MONX11151



Mouse anti-Cytokeratin 8, clone TS1 (monoclonal) Clone no. TS1

Additional info

Analyte Specific Reagent. Analytical and performance characteristics are not established.

References	1.	Johansson A et al. Cancer Research. 59: 48–51 (1999
	2	Sundström BE et al. The Journal of Histochemistry and Cytochemistry. 37 (12):
	3.	-
	4.	-
	5.	-

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021