

Rabbit anti-p53 Protein (CM5) (polyclonal)

Clone no. - MONXtra

Product name Rabbit anti-p53 Protein (CM5) (polyclonal)

Host Rabbit

Applications IHC-P (1:200), WB (1:2000-1:4000)

Species reactivity human

Conjugate -

Immunogen Recombinant mouse p53 protein

Isotype -

Clonality Polyclonal

Clone number -

Size 0.5ml

Concentration n/a

Format -

Storage buffer PBS with 1% BSA and sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MONX11159



Rabbit anti-p53 Protein (CM5) (polyclonal)

Clone no. - MONXtra

Additional info

The accumulation of p53 protein in response to genotoxic stress in vitro is well established and appears to induce growth arrest and apoptosis by the transcriptional regulation of other genes and possibly by other direct mechanisms.

References 1. Botchkarev VA et al. American Journal of Pathology. 158 (6): 1913–1919 (2001)

2 Midgley CA et al. Journal of Cell Science. 108: 1843–1848 (1995)

3. Milner JO et al. Oncogene. 8: 2001–2008 (1993

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES