Product datasheet

PS006



MONOSAN

Polyclonal Rabbit anti-rat IFN-gamma Clone no. -

Product name	Polyclonal Rabbit anti-rat IFN-gamma	
Host	Rabbit	
Applications	IHC-fr,FUNC,ELISA,WB	
Species reactivity	rat	
Conjugate	-	
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers	
lsotype	lg	
Clonality	Polyclonal	
Clone number	-	
Size	1 ml	
Concentration	100 ug/ ml	
Format	-	
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS006

Polyclonal Rabbit anti-rat IFN-gamma Clone no. -

MONOSAN[°]

MONOSAN

Additional info

The antibody reacts with rat Interferon gamma (IFN-gamma) of both natural and recombinant origin. IFN-gamma is a pluripotent cytokine with important pro-inflammatory functions. The antibody inhibits the biological activity of natural and recombinant rat IFN-gamma. The antigen specificity was further assessed by ELISA and Western blotting. No cross-reactivities with other cytokines have been detected.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com