Product datasheet

PS014



MONOSAN

Polyclonal Rabbit anti-Mouse Interleukin-6 Clone no. -

Product name	Polyclonal Rabbit anti-Mouse Interleukin-6	
Host	Rabbit	
Applications	FUNC,ELISA,WB	
Species reactivity	mouse, rat	
Conjugate	-	
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers	
lsotype	lg	
Clonality	Polyclonal	
Clone number	-	
Size	1 ml	
Concentration	100 ug/ ml	
Format	-	
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

Polyclonal Rabbit anti-Mouse Interleukin-6 Clone no. -

MONOSAN[°]

MONOSAN

Additional info

The antibody reacts with natural and recombinant mouse interleukin 6 (IL-6) as assessed by ELISA. The antibody inhibits the biological activity of natural and recombinant mouse IL-6 as determined with the B9 cell bio-assay. The antibody cross-reacts with natural rat IL-6. IL-6 is a pluripotent 20-22 kDa cytokine which plays a role in the pathophysiology of severe infection and regulates the immune response, acute phase reaction and hematopoiesis. IL-6 plays a critical role in B-cell differentiation to plasma cells and is a potent growth factor for plasmacytoma and myeloma. Continuous IL-6 gene expression makes an essential contribution to a multistep oncogenesis of plasma cell neoplasia. IL-6 is a very useful culture supplement for the generation of a high number of antibody-producing hybridomas.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com