Product datasheet PS016



## Polyclonal Rabbit anti-Interleukin-10

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-Interleukin-10

**Host** Rabbit

**Applications** IHC-fr,FUNC,ELISA,IF,IHC-P,WB

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lg

**Clonality** Polyclonal

Clone number -

Size 1 ml

Concentration 100 ug/ ml

Format -

Storage buffer PBS with 0.1% BSA and 0.02% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS016



Polyclonal Rabbit anti-Interleukin-10

Clone no. - MONOSAN

## Additional info

The antibody reacts with human Interleukin-10 (IL-10) of both natural and recombinant origin. IL-10 is a pluripotent cytokine with important immunosuppressive actions: it can inhibit cytokines involved in the Th1 response such as IL-2, Interferon gamma and also inhibits the production of pro-inflammatory cytokines such as IL-1, IL-6, IL-8 and TNF-alpha. The antibody inhibits the biological activity of natural and recombinant human IL-10. The antigen specificity was further assessed by ELISA. No cross reactivities with other cytokines have been detected.

References 1. -

5.

2 -

3. -

4. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES