

Product datasheet PS018

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-ICAM-1, CD54 (Biotinylated)

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-ICAM-1, CD54 (Biotinylated)
Host	Rabbit
Applications	IHC-fr,ELISA,IP,WB
Species reactivity	human, baboon, chimpanzee, cynomolgus monkey, rhesus monkey
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	Ig
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	50 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-ICAM-1, CD54 (Biotinylated)

Clone no. -

MONOSAN

**Additional info**

The antibody reacts with both intact human ICAM-1, CD54 and with soluble human ICAM-1. The antibody reacts also with ICAM-1 present on chimpanzee, rhesus monkey, cynomolgus monkey and baboon tissues.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES