

Product datasheet PS019

MONOSAN[®]

Polyclonal Rabbit anti-Interferon-Gamma

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Interferon-Gamma
Host	Rabbit
Applications	IHC-fr,FUNC,ELISA,WB
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	Ig
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Interferon-Gamma

Clone no. -

MONOSAN

Additional info

The antibody reacts with human gamma IFN (IFN-gamma) of both natural and recombinant origin. IFN-gamma is a pluripotent cytokine with important pro-inflammatory functions. The antibody inhibits the biological activity of natural and recombinant human IFN-gamma. The antigen specificity was further assessed by ELISA and Western blot. No cross reactivities with other cytokines have been detected.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES