

Product datasheet PS033

MONOSAN[®]

Polyclonal Rabbit anti-Keratin

Clone no. -

MONOSAN

Product name Polyclonal Rabbit anti-Keratin

Host Rabbit

Applications IHC-fr, IHC-P (1:50)

Species reactivity Human, Mouse

Conjugate -

Immunogen Human skin

Isotype IgG

Clonality Polyclonal

Clone number -

Size 250 µl

Concentration n/a

Format -

Storage buffer Contains 0.01% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021

Revision date
22-12-2021

Polyclonal Rabbit anti-Keratin

Clone no. -

MONOSAN

Additional info

The antiserum reacts positively with all stratified squamous epithelia and epidermal appendages such as hairfollicles, sebaceous glands. Also positive on epithelia of the urinary tract, stomach, intestine, uterine endometrium and prostate. Identifies mesotheliomas, squamous cell carcinomas, adenocarcinomas, transitional cell carcinomas and anaplastic carcinomas.

Recommended for positive control: Squamous cell carcinoma, skin, cervix.

References

1. Ramaekers F et al. A Pathol Anat Histopathol 1985;408:127-142
2. Yanagi T et al. Cancer Res 2018;78:6795-6806
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES