

Product datasheet PS038

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Interleukin-8

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Interleukin-8
Host	Rabbit
Applications	FUNC,ELISA,IP,WB
Species reactivity	human, cynomolgus monkey, rhesus monkey
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	Ig
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Interleukin-8

Clone no. -

MONOSAN

**Additional info**

The polyclonal antibody reacts with human natural and recombinant Interleukin-8 (IL-8) as assessed by ELISA. The antibody inhibits the biological activity of human native and recombinant IL-8. The antibody cross reacts with rhesus and cynomolgus natural IL-8.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES