

Product datasheet PS052

MONOSAN[®]

Polyclonal Rabbit anti-Mouse TNF-alpha

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Mouse TNF-alpha
Host	Rabbit
Applications	IHC-fr,FC,FUNC,ELISA,IHC-P,WB
Species reactivity	mouse, rat
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	Ig
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Mouse TNF-alpha

Clone no. -

MONOSAN

Additional info

Tumor Necrosis Factor alpha (TNF-alpha) is a cytokine which has diverse immunomodulatory, anti-tumor and toxic effects. TNF-alpha has been detected in diverse inflammatory status and appears to be a critical mediator in the lethality of septic shock. Furthermore, TNF-alpha has also been found in inflammatory foci such as synovial effusions in rheumatoid arthritis, systemic circulation in septic shock, parasitemia and rejection of renal transplants. The antibody reacts with both rat and mouse natural and recombinant TNF-alpha and recognizes membrane and receptor bound TNF-alpha. The antibody shows neutralizing activity.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES