

Product datasheet PS053

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-CD120a, TNF-RI (p55/p60)

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-CD120a, TNF-RI (p55/p60)
Host	Rabbit
Applications	FUNC,ELISA,IP,WB
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	Ig
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-CD120a, TNF-RI (p55/p60)

Clone no. -

MONOSAN

**Additional info**

The antibody reacts with the extra-cellular part of the TNF-RI and with the soluble receptor. TNF-RI is present on most cell types and is considered to play a prominent role in cell stimulation by TNF-alpha. Induction of cytotoxicity and other functions are mediated largely via TNF-RI.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES