

Polyclonal Rabbit anti-Rat Nitric oxide synthase

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Rat Nitric oxide synthase
Host	Rabbit
Applications	IHC-fr (1:1,000-1:5,000)
Species reactivity	Cat, Human, Rat
Conjugate	-
Immunogen	Synthetic peptide from the C-terminus of rat cerebellar NO-synthase, coupled to BSA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Rat Nitric oxide synthase

Clone no. -

MONOSAN

Additional info

Nitric oxide synthase is an enzyme catalizing the synthesis of NO from L-arginine. The antiserum was raised against a synthetic peptide of rat cerebellar NOS, that shows no homology with other related proteins.

Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining.

Positive control: frozen sections of rat colon.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES