

Product datasheet PS080

MONOSAN[®]

Polyclonal Rabbit anti-ACTH

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-ACTH
Host	Rabbit
Applications	(IHC-fr, IHC-P)1:100 – 1:1200
Species reactivity	Human, Pig
Conjugate	-
Immunogen	Synthetic ACTH
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	250 µl
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-ACTH

Clone no. -

MONOSAN

Additional info

ACTH may occur in pituitary adenomas, medullary thyroid carcinomas, bronchial carcinoids, pancreatic islet cell tumors and gastric carcinoids.

Antigen localization: cytoplasm.

Recommended for positive control: pituitary gland

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES