

Product datasheet PS085

MONOSAN[®]

Polyclonal Rabbit anti-Gastrin

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Gastrin
Host	Rabbit
Applications	IHC-fr, IHC-P (1:50-1:100), IF
Species reactivity	Human, Mammals
Conjugate	-
Immunogen	Human gastrin 17 conjugated to BSA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	250 µl
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Gastrin

Clone no. -

MONOSAN

Additional info

Gastrin-secreting cells are numerous in the antrum and a few are found in the proximal duodenum. The antibody can be used for the diagnosis of gastrin-producing tumors which are mainly found in the pancreas and occasionally in the stomach and the duodenum.
Positive control: Antrum (antigen localization: cytoplasmic, extracellular)

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES