

Product datasheet PS088

**MONOSAN**<sup>®</sup>

Polyclonal Rabbit anti-Insulin  
Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Insulin
Host	Rabbit
Applications	IHC-fr, IHC-P (1:50-1:100)
Species reactivity	Human, Rat
Conjugate	-
Immunogen	Porcine insulin
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	250 µl
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Insulin

Clone no. -

MONOSAN

**Additional info**

Insulin is produced by the B-cells of the pancreatic islets. The antibody can be used for the diagnosis of pancreatic B-cell tumors. <br />Antigen localization: cyto-plasm, extracellular.<br />Positive control: Pancreas.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES