Product datasheet PS089



Polyclonal Rabbit anti-Bombesin

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-Bombesin

Host Rabbit

Applications IHC-fr, IHC-P (1:50)

Species reactivity Frog, Human, Mouse, Rat

Conjugate -

Immunogen Synthetic frog bombesin conjugated to bovine thyroglobulin

Isotype -

Clonality Polyclonal

Clone number -

Size 250 µl

Concentration n/a

Format -

Storage buffer Contains 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS089



Polyclonal Rabbit anti-Bombesin

Clone no. - MONOSAN

Additional info

This polyclonal antibody can be used for the study of endocrine differentiation in tumors of the respiratory tract.

Positive control: Gut (peripheral nerves), fetal lung

References 1. -

2 -

3. -

4. -5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES