

Product datasheet PS101P

MONOSAN[®]

Polyclonal Goat anti-Human IgM, HRP

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Human IgM, HRP
Host	Goat
Applications	WB, ELISA, IHC
Species reactivity	human
Conjugate	HRP
Immunogen	IgM
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.01 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS101P

MONOSAN[®]

Polyclonal Goat anti-Human IgM, HRP

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/ml in the buffer (0.1 M sodium chloride, 0.01M sodium phosphate pH 7.5, 50% glycerol, 5 mg/ml BSA). Specific to human IgM (heavy chain). No cross-reactivity towards human IgG, IgA and other human plasma proteins

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES