Product datasheet PS203



Polyclonal Rabbit anti-alpha-MSH

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-alpha-MSH

Host Rabbit

Applications IHC-fr, IHC-P (1:200-1:400)

Species reactivity Pig, Rat

Conjugate -

Immunogen Synthetic alpha-MSH (Sigma Chem. Co., USA) conjugated to BSA

Isotype -

Clonality Polyclonal

Clone number -

Size 50 µl (lyoph.)

Concentration n/a

Format -

Storage buffer Lyophilized; reconstitute in 100 µl dist. water

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS203



Polyclonal Rabbit anti-alpha-MSH

Clone no. - MONOSAN

Additional info

Specific for alpha-melanocyte stimulating hormone. Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining.

strong>alpha-MSH is found in corticotrophs of the anterior pituitary and may also occur in brain neurons.

br>Positive control: Bouinfixed paraffin sections of pig pituitary.

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES