

Polyclonal Rabbit anti-Neuromedin U-8

Clone no. -

MONOSAN

---

Product name	Polyclonal Rabbit anti-Neuromedin U-8
Host	Rabbit
Applications	IHC-fr, IHC-P (1:1,000-1:1,500)
Species reactivity	Chicken
Conjugate	-
Immunogen	Synthetic neuromedin U-8, porcine sequence, conjugated to BSA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Neuromedin U-8

Clone no. -

MONOSAN

---

**Additional info**

Neuromedin U was found in porcine spinal cord. Neuromedin U-8 is contained within a larger polypeptide form neuromedin U-25. Neuromedin U is present in central and peripheral neurons, particularly in the enteric nervous system. <br>Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining.

Positive control: frozen sections of chicken small intestine.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES