

Product datasheet PS208

MONOSAN[®]

Polyclonal Rabbit anti-Motilin

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Motilin
Host	Rabbit
Applications	IHC-fr, IHC-P (1:500-1:5,000)
Species reactivity	Cat, Human, Pig
Conjugate	-
Immunogen	Pure porcine motilin (Prof. V. Mutt), conjugated to BSA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Motilin

Clone no. -

MONOSAN

Additional info

Specificity: absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining.

Positive control: formalin-fixed paraffin sections of pig duodenum.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES