Product datasheet

PS211



MONOSAN

Polyclonal Rabbit anti-VAChT Clone no. -

Product name	Polyclonal Rabbit anti-VAChT
Host	Rabbit
Applications	IHC-fr, IHC-P (1:800-1:1,200)
Species reactivity	Human, Rat
Conjugate	_
Immunogen	Synthetic peptide from the C-terminus of rat VAChT conjugated to BSA
lsotype	-
Clonality	Polyclonal
Clone number	_
Size	50 μl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet



Polyclonal Rabbit anti-VAChT Clone no. -

MONOSAN

Additional info

The antiserum against the vesicular acetylcholine transporter is a unique immunohistochemical marker for cholinergic nerves, more specific than the commonly used acetylcholine esterase (AchE), since it does not react with postsynaptic neurons, and is more sensitive than choline acetyltransferase (ChAT). The antiserum recognizes VAChT both in the CNS and PNS of rat and mouse.

br>Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes staining.

Positive control: frozen sections of rat small intestine.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com