

Product datasheet PS217

MONOSAN[®]

Polyclonal Rabbit anti-Secretin

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Secretin
Host	Rabbit
Applications	IHC-fr, IHC-P (1:400-1:800)
Species reactivity	Human, Pig
Conjugate	-
Immunogen	Pure porcine secretin (Prof. V. Mutt)
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Secretin

Clone no. -

MONOSAN

Additional info

Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining, while glucagon, GIP, PP and VIP do not.
Positive control: frozen sections of rat duodenum.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES