

Product datasheet PS224

MONOSAN[®]

Polyclonal Rabbit anti-Neurotensin

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Neurotensin
Host	Rabbit
Applications	IHC-fr, IHC-P (1:1,000-1:3,000)
Species reactivity	Cat
Conjugate	-
Immunogen	Synthetic bovine neurotensin (Beckman, Ca. USA), conjugated to BSA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Neurotensin

Clone no. -

MONOSAN

Additional info

Absorption with 10-100 ug neurotensin per ml diluted antiserum abolishes the staining while glucagon, insulin, gastrin, GIP and VIP do not.

Positive control: Bouin-fixed paraffin or Stefanini-fixed frozen sections of cat ileum.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES