

Product datasheet PS229

MONOSAN[®]

Polyclonal Rabbit anti-PYY

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-PYY
Host	Rabbit
Applications	IHC-fr, IHC-P (1:1,000-1:5,000)
Species reactivity	Human, Rat
Conjugate	-
Immunogen	Synthetic porcine peptide YY (PYY) (Peninsula, code 7173)
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-PYY

Clone no. -

MONOSAN

Additional info

Specificity: absorption with 10-100 ug PYY per ml diluted antiserum abolishes the staining, while NPY and PP do not.

Positive control: Bouin-fixed paraffin sections of rat colon and frozen sections of rat colon.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES