

Product datasheet PS230

MONOSAN[®]

Polyclonal Rabbit anti-Prepro-NPY 68-97 (C-PON)

Clone no. -

MONOSAN

Product name Polyclonal Rabbit anti-Prepro-NPY 68-97 (C-PON)

Host Rabbit

Applications IHC-fr (1:400-1:800)

Species reactivity Human, Rat

Conjugate -

Immunogen Synthetic Prepro-NPY 68-97 (C-PON)

Isotype -

Clonality Polyclonal

Clone number -

Size 50 µl (lyoph.)

Concentration n/a

Format -

Storage buffer Lyophilized; reconstitute in 100 µl dist. water

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021

Revision date
22-12-2021

Polyclonal Rabbit anti-Prepro-NPY 68-97 (C-PON)

Clone no. -

MONOSAN

Additional info

Neuropeptide Y is a peptide belonging to the PP-family and occurs in neurons and adrenal medullary cells. Antisera raised against NPY often cross-react with PP and PYY. This antiserum was raised using a synthetic peptide from the prepro-NPY sequence which has no homologies to PP and NPY.

Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining.

Positive control: Stefanini-fixed sections of rat small intestine.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES