Product datasheet PS241



Polyclonal Rabbit anti-Rat Pancreastatin

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-Rat Pancreastatin

Host Rabbit

Applications IHC-fr, IHC-P (1:200-1:1,000)

Species reactivity Rat

Conjugate -

Immunogen Synthetic chromogranin A (aa 306-314), rat

Isotype -

Clonality Polyclonal

Clone number -

Size 50 μl (lyoph.)

Concentration n/a

Format -

Storage buffer Lyophilized; reconstitute in 100 µl dist. water

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS241



Polyclonal Rabbit anti-Rat Pancreastatin

Clone no. - MONOSAN

Additional info

Pancreastatin is a fragment of chromogranin A and is produced by proteolytic processing of chromogranin A in several peptide hormone-producing cells, such as pancreatic islet cells and gut endocrine cells, as well as tumors arising from these cells.
br> Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes the staining. Positive control: Bouin-fixed paraffin sections of rat pancreas.

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES